



Python and GUI Implementation for Internal Optimization and Increased Usability of *in vivo* Glucose Sensing System



Kathryn DiPippo^{1,2}, Alex Werth¹, Michelle Zhang^{1,3}, Phillip Heucke^{1,4}, Jessica Doyle^{1,5}, Claire Gmachl¹

1 – Princeton University Department of Electrical Engineering, Princeton, NJ, 08540, USA

2 – Permanent Address: Rensselaer Polytechnic Institute, Troy, NY, 12180, USA





3 – Permanent Address: Smith College, Northampton, MA, 01063, USA

Industry / Practitioner Collaborators: Daylight Solutions Inc. Zurich Instruments

4 – Permanent Address: Bucknell University, Lewisburg, PA, 17837, USA

5 – Permanent Address: Hunterdon Central Regional High School, Flemington, NJ, 08822, USA

email: dipippo@princeton.edu / dipipk@rpi.edu



- GUI started up and populated with data
- Number of spectral runs are taken
- Various statistical analyses applied:
 - Wavelet transform
 - Savitzky-Golay Filtering with Differentiation
 - Machine Learning Prediction Algorithm
- Final predicted glucose concentration obtained



GUI – CLINICAL AND SOLUTIONS



File Edit

1. Store Data:

Select the Folder to Store Data:

2. Signal Processing:

Number of Runs:

START Frequency:

STOP Frequency:

Number of Scans:

3. Basic Information:

Name:

Concentration: mg/L

Date:


4. Real-Time Calculation:

☐ Remove Outlier?

Final Number of Runs:



























☐ Display Plot?

5. Run:









ge, Northampton, MA, 01063, USA
 nic Institute, Troy, NY, 12180, USA
 nt Solutions Inc. Zurich Instruments

-  _DS_Store
-  032615.mat
-  addManually.m
-  addPatient.m
-  clinicGUI.fig
-  clinicGUI.m
-  clinicGUIUpdated.fig
-  combineSets.m
-  findminstd2.m
-  focusing_lens_test_052115.jpg
-  FrozenComputer.m
-  getDataStruct.m
-  getOtherDataStruct.m
-  InitLaser.m
-  LaserPlusLIF.m
-  LaserPlusLIF2.m
-  lens_2_052015.mat
-  lens_52015.mat
-  LoadMat.m
-  matlab.mat
-  pfile.m
-  pfile1.m
-  pfile2.m
-  test_function.m
-  Thumbs.db
-  Zurich_async_SINGLE.m



-  clinical.py
-  combinedGUI.py
-  GUI_BACKGROUND.gif
-  solutions.py

MatLab files needed prior to conversion

